

# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Version Number: 9 Issue Date:16 April 2019

## SECTION 1: Identification of the substance/mixture and company/undertaking

1.1 **Product Identifier** 

> **Product Name** Clinell Universal Wipes

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified Use Wipes for surface disinfection and cleaning of non-invasive

medical devices.

1.3 Details of the supplier of the safety data sheet

Supplier GAMA Healthcare Ltd

2 Regal Way Watford WD24 4YJ United Kingdom

Tel: +44 (0) 207 993 0030

Email: info@gamahealthcare.com

**Emergency telephone number** 

Tel: +44 (0) 207 9930 035

#### **SECTION 2: Hazards identification**

#### 2.1 Classification of the substance or mixture

Classification according Mixture not classified as hazardous to Regulation (EC) No

1272/2008

#### 2.2 Label Elements

Contains PHMB. May produce an allergic reaction. Wear appropriate PPE

#### 2.3 Other hazards

Not applicable

#### **SECTION 3: Composition/information on ingredients**

#### 3.2 Mixtures

Declarable components	Conc. (%)	EC No.	CAS No.	Classification of individual components under Regulation EC No1272/2008	
Benzalkonium chloride	≤0.5	270-325-2	68424-85-1	Skin Corr 1B (H314) Acute Tox 4 (H302, H312) Aguatic Acute 1 (H400)	
Didecyl dimethyl ammonium chloride	≤0.5	230-525-2	7173-51-5	Acute Tox 4 (H302) Skin Corr 1B (H314)	
Polyhexamethylene biguanide (PHMB)	≤0.10	NA	27083-27-8	Acute tox 4 (H302) Skin sens 1B (H317) Eye dam 1 (H318) Carc. 2 (H351) STOT RE 1 (H372)	



# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Issue Date:16 April 2019 Version Number: 9

				Aquatic acute 1 (H400) Aquatic chronic (H410)
--	--	--	--	--

Other components:

Water >75 Additives Each <1

#### **SECTION 4: First aid measures**

#### 4.1 Description of first aid measures

#### Inhalation

Acute effects following exposure to this product via the inhalation route are not anticipated during normal handling and use.

#### Skin

This product is not intended for skin use, wear PPE. If irritation develops, seek medical advice.

Although this product contains components classified as corrosive and sensitising to skin, due to the high volume of water also present in the formulation, the dilution effect means the classification of the formulation through CLP does not result in the hazard being carried through to the product.

A toxicological risk assessment considers this product unlikely to causes significant dermal irritation, sensitisation or delayed hypersensitivity.

# Eye

The product contains components classified as damaging to eyes. Due to the high volume of water also present in the formulation, the dilution effect means that the classification through CLP does not result in the hazard being carried through to the product.

Nevertheless, should eye irritation be experienced, this effect would likely be transient. But should symptoms persist, seek medical advice.

#### Ingestion

This product is for external use only and should be kept away from children. No adverse effects are anticipated from the formulation via the oral route during normal handling and use of the product.

# **4.2** Most important symptoms and effects, both acute and delayed This product contains PHMB which may cause an allergic reaction.

# **4.3** Indication of any immediate medical attention and special treatment needed Treat symptoms as they occur.



# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Issue Date:16 April 2019 Version Number: 9

#### **SECTION 5: Firefighting measures**

#### 5.1 Extinguishing media

Water spray, carbon dioxide, dry chemical and foam are compatible with the product. No unsuitable extinguishing media are known.

#### 5.2 Special hazards arising from the substance of mixture

The product is water based, therefore not flammable or explosive.

## 5.3 Advice for fire fighters

Fire fighters should wear an approved self-contained breathing apparatus and full protective clothing.

#### **SECTION 6: Accidental release measures**

# **6.1** Personal precautions, protective equipment and emergency procedures None anticipated or expected to be required.

#### 6.2 Environmental precautions

None anticipated or expected to be required.

#### 6.3 Methods and material for containment and cleaning up

None anticipated or expected to be required.

#### 6.4 Reference to other sections

For recommended personal protective equipment see Section 8.

#### **SECTION 7: Handling and storage**

#### 7.1 Precautions for safe handling

For prolonged use wear gloves to avoid drying of the skin.

#### 7.2 Conditions for safe storage, including any incompatibilities

Store in a cool, dry, well-ventilated place, away from direct sunlight. Do not allow to freeze. Keep container closed when not in use.

#### 7.3 Specific end use

See directions for use on pack.

Identified in Section 1.2



# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Issue Date:16 April 2019 Version Number: 9

#### **SECTION 8: Exposure controls/personal protection**

#### 8.1 Control Parameters

**EU Limit:** 

No applicable EU occupational exposure limit values

## 8.2 Exposure controls

Engineering controls

None anticipated or expected to be required.

Personal protective equipment

For prolonged use, wear gloves.

Environmental exposure controls

None anticipated or expected to be required.

#### **SECTION 9: Physical and chemical properties**

#### 9.1 Information on basic physical and chemical properties

Appearance Moist non woven wipe

Odour Slight green tea perfume

Odour threshold Not available

pH 5-8

Melting/freezing point Ca. 0°C Initial boiling point/range Ca. 100°C

Flash point Not determined: water based product

Evaporation rate Not determined: water based product

Flammability (solid, gas) Not determined: water based product

Flammability or explosive limits 
Not determined: water based product



# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Issue Date:16 April 2019 Version Number: 9

Vapour pressure 24 mmHg (25°C) (water)

Relative density Not determined: water based product

Solubility Not determined: water based product

Partition coef No data available

Auto-ignition temperature Not determined: water based product

Decomposition temperature Not determined: water based product

Viscosity Not determined: water based product

Explosive properties Not determined: water based product

Oxidising properties Not determined: water based product

**9.2 Other information** Not available

#### **SECTION 10: Stability and reactivity**

#### 10.1 Reactivity

Contact with ionic substances for example oils and dyes, may reduce effectiveness of the product. Contact with oxidising agents should be avoided.

#### 10.2 Chemical stability

The product is considered stable under normal ambient storage and handling conditions or temperature and pressure.

#### 10.3 Possibility of hazardous reactions

No hazardous reactions anticipated

#### 10.4 Conditions to avoid

None known

#### 10.5 Incompatible materials

Oxidizing agents and anionic formulations.

#### 10.6 Hazardous decomposition products

None known.

#### **SECTION 11: Toxicological information**

This preparation has undergone toxicology risk assessment.



# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Issue Date:16 April 2019 Version Number: 9

# 11.1 Information of toxicological effects

Acute toxicity

Not likely to be acutely toxic.

Irritancy

Not likely to cause significant dermal irritation.

Corrosivity

No risk of dermal corrosivity identified under normal handling and use.

Sensitisation

Not likely to cause significant sensitisation or delayed hypersensitivity.

Repeated dose toxicity

No data available on the repeat dose toxicity of this product.

Carcinogenicity

No data available on the carcinogenicity of this product.

Mutagenicity

None of the components have exhibited confirmed mutagenic characteristics in the evaluation of their toxicity to date.

Toxicity for reproduction

Not data available on toxicity for reproduction of this product.

#### **SECTION 12: Ecological information**

Ecotoxicological data has not been determined specifically for this product. Based on classification of the formulation through CLP, the environmental hazards are not carried through to the product.

#### **12.1 Toxicity**

Components are classified as toxic to the environment but are not present in the formulation at sufficient levels. The hazard is not carried through to the product.

# 12.2 Persistence and degradability

Two components of the formulation (DDAC and BAC) have been found to readily biodegrade in OECD 301D closed bottle tests. However, PHMB was found not to be readily biodegradable under the same protocol.

## 12.3 Bioaccumulative potential



# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Issue Date:16 April 2019 Version Number: 9

Due to the distribution coefficient of n-octonal/water, accumulation in organisms is not expected.

#### 12.4 Mobility soil

No information available on mobility of active substance in soil.

#### 12.5 Results of PBT and vPVB assessment

The formulation does not contain substances that meet the PBT or vPvB criteria of REACH annex XIII.

#### 12.6 Other adverse effects

No additional information available

#### **SECTION 13: Disposal considerations**

#### 13.1 Waste treatment methods

This product may be disposed of landfill, or by incineration. Disposal must be in accordance with current national and local regulations.

In the Healthcare Industry, chemical residues, biocides and infectious substances generated as a result of medical and nursing care may require classification as hazardous waste.

Waste disposal is regulated in the EC member countries through corresponding laws and regulation. In the UK, we recommend that you consult the List of Wastes available through the Environment Agency. In other countries, contact either the authorities or approved waste disposal companies for advice on disposal of used waste.

#### **SECTION 14: Transport Information**

Not classified for transport

#### **SECTION 15: Regulatory information**

## 15.1 Safety, health and environmental regulations/legislation specific for the mixture

This product is classified under the Classification, Labelling and Packaging of Substances and Mixtures (EC) No 1272/ 2008 it contains substances which are notified and under the Biocidal Products Regulation (EU). No 528/2012.

#### 15.2 Chemical safety assessment

Not applicable

#### **SECTION 16: Other Information**



# **Clinell Universal Wipes**

According to Regulation (EU) No 453/2010

Issue Date:16 April 2019 Version Number: 9

Revisions

Currently in its eighth revision.

Basis of classification

The mixture is self-classified on the basis of available information on the ingredients.

This safety data sheet was compiled using the ECHA Guidance on the compilation of Safety Data Sheets, Version 1.1 December 2011.

Disclaimer: This information is furnished without warranty express or implied, except that it is accurate to the best of our knowledge.

It relates to the specific material designated herein, and does not relate to the use in combination with any other material or in any process. We assume no legal responsibility for use or reliance upon this information.

# **cline**Universal Sanitising Wipes Efficacy Data

Effective Against Viruse	Test	Kill Time
COVID-19 (SARS-CoV-2)	EN14476	30 sec
nfluenza H1N1	EN14476	2 min
nfluenza H3N2	EN14476	2 min
Hepatitis B	ASTM E1052	1 min
Hepatitis C	EN14476	1 min
Herpes Simplex Type 1	ASTM E1053	2 min
HIV	EN14476	30 sec
MERS-CoV	EN14476	1 min
Measles vius (Canine distemper virus)	EN14476	1 min
Mumps virus (Rubulavirus)	EN14476	1 min
Murine coronavirus	EN14476	30 sec
Murine Hepatitis Virus (Betacoronavirus surrogate)	EN14476	30 sec
Norovirus	EN14476	1 min
Rotavirus	EN14476	1 min
Vaccinia virus	EN14476	15 sec
Bacte	ria	
	EN13727	10 sec
Acinetobacter baumannii	EN13727	1 min
	EN16615	2 min
Burkholderia cepacia	EN13727	10 sec
Enterococcus faecalis	EN13727	10 sec
	EN16615	2min
Enterococcus faecium	EN13727	10 sec
	EN13697	2 min
Enterococcus hirae	EN13727	10 sec
	EN16615	10 sec
Escherichia coli	EN13697	2 min
ESCHERCING COIL	EN16615	2 min
Escherichia coli k12	EN13727	10 sec
Klebsiella pneumoniae	EN13727	10 sec
Klebsiella pneumoniae (CPE)	EN13727	10 sec
Klebsiella pneumoniae (ESBL)	EN16615	2 min
Legionella pneumophila	EN13727	1 min
	EN13727	10 sec
Methicillin-resistant Staphylococcus aureus	EN13727	1 min
	EN16615	2 min
	EN13697	2 min
Pseudomonas aeruginosa	EN13727	10 sec
	EN16615	10 sec
Salmonella typhimurium	EN13727	1 min
Serratia marsescens	EN13727	10 sec
	EN13697	5 min
Staphylococcus aureus	EN13727	10 sec
	EN166515	10 sec
Vancomycin resistant enterococcus faecium (VRE)	EN13727	10 sec
Vancomycin resistant enterococci (VRE)	EN13727	1 min
<u> </u>	EN16615	2 min
Fung		lie :
Aspergillus brasiliensis	EN13697	15 min
Mycoba:		I o main
Mycobacterium var bovis  Yea:	EN14348	2 min
Teas	EN13624	1 min
Candida albicans	EN13624 EN13624	2 min
Canalaa albicans	EN 13624 EN 16615	
Candida auris	EN13624	10 sec
Candida auris Candida auris Japanese clade	EN13624	1 min
Candida auris South African clade	EN13624	1 min
Candida auris South Airican clade  Candida auris South Korean clade	EN13624 EN13624	1 min
Veterinary o		[ T 111111
Bordetella bronchiseptica	EN1656	1 min
Candida albicans	EN1657	1 min
Enterococcus hirae	EN1656	1 min
	EN1656	1 min
entospira interrogans	I LINIUJO	[ 1 1 1 m 1
<u> </u>		1 min
Methicillin-Resistant Staphylococcus pseudintermedius	EN1656	1 min
Leptospira interrogans  Methicillin-Resistant Staphylococcus pseudintermedius  Proteus vulgaris  Pseudomonas geruginosa	EN1656 EN1656	1 min
Methicillin-Resistant Staphylococcus pseudintermedius	EN1656	